

PAJUNK®



Pain Management Product Catalogue

Table of Contents

1 | INJECTIONS

<i>SonoMSK</i>	4
<i>RGN</i>	5
<i>TuohySono</i>	6
<i>Tuohy II</i>	7
<i>LOR</i>	8

2 | RF ABLATION

<i>RFTL</i>	12
<i>RFTL Curved</i>	13
<i>RFTL Sono</i>	14
<i>PulsedRF</i>	15
<i>AlliClip</i>	16

Injections


*Ultrasound and Fluoroscopy Guided Steroid
& Anaesthetic Injection Needles*



SonoMSK

Sharp, Echogenic MSK and Chronic Pain Injection Needle

- 1 Echogenic Needle Tip
- 2 Cornerstone Reflectors
- 3 Sharp Back Cut Bevel
- 4 Depth Markings
- 5 Pre-attached & Removable Injection Tubing



Product advantages at a glance:


- Reliable and optimised needle echogenicity at higher angle
- Optimised ultrasound visibility of needle shaft
- Reduced pain during puncture and needle placement
- Smooth gliding properties
- Easier and more precise needle maneuverability
- Convenience

Description	Item no.	PU
27G x 37 mm (1 1/3")	1191-4B037	10
27G x 50 mm (2")	1191-4B050	10
25G x 37 mm (1 1/3")	1191-4C037	10
25G x 50 mm (2")	1191-4C050	10
25G x 70 mm (2 3/4")	1191-4C070	10
25G x 90 mm (3 1/2")	1191-4C090	10
22G x 25 mm (1")	1191-4E025	10
22G x 40 mm (1 5/8")	1191-4E040	10
22G x 70 mm (2 3/4")	1191-4E070	10
22G x 90 mm (3 1/2")	1191-4E090	10
22G x 100 mm (4")	1191-4E100	10
22G x 120 mm (4 3/4")	1191-4E120	10
22G x 150 mm (6")	1191-4E150	10

RGN

Hubless Contrast Injection Needle

- 1 Hubless Design
- 2 Fixed Injection Tubing
- 3 Sharp Facet Tip
- 4 Smooth Needle Surface



Product advantages at a glance:


- Reduced pain during puncture and needle placement
- Smooth gliding properties
- Reduced shadowing and improved visibility of target area
- Convenience due to closed system

Description	Item no.	PU
23G x 60 mm (2 1/4")	001160-10	10
22G x 100 mm (4")	001160-20	10
21G x 150 mm (6")	001160-40	10

TuohySono

Echogenic, Non-Stimulating Single Shot Fascial Plane Block Needle

- 1 Echogenic Needle Tip
- 2 Cornerstone Reflectors
- 3 Depth Markings
- 4 Precision-Ground Metal Stylet
- 5 Detachable Plastic Wings



Product advantages at a glance:

- Reliable and optimised needle echogenicity at higher angle
- Optimised ultrasound visibility of needle shaft
- Minimised tissue coring
- Smooth gliding properties
- Optimum handling for different gripping techniques


Description	Item no.	PU
21G x 90mm	1085-4F090	25
• 18G x 50 mm (2")	1085-4K050	25
18G x 90 mm (3 1/2")	1085-4K090	25
18G x 120 mm (4 3/4")	1085-4K120	10
17G x 90 mm (3 1/2")	1085-4M090	25

- Small patient product

Tuohy II

Non-Stimulating Single Shot Fascial Plane Block Needle

- 1 Depth Markings
- 2 Precision-ground Plastic Stylet
- 3 Integrated Ergonomic Wings



Product advantages at a glance:

- Minimised tissue coring
- Smooth gliding properties
- Optimum handling for different gripping techniques

Description	Item no.	NRFit® Item no.	PU
19G x 90mm (3 1/2")	1170-6I091	1176-6I091	10
18G x 80mm (3 1/4")	1170-6K081	1176-6K081	10
18G x 90mm (3 1/2")	1170-6K091	1176-6K091	10
18G x 120mm (4 3/4")	1170-6K120	1176-6K120	10
17G x 90mm (3 1/2")	1170-6M091	1176-6M091	10
17G x 120mm (4 3/4")	1170-6M120	1176-6M120	10

LOR Syringe



Description	Item no.	PU
Plastic LOR Syringe 10 ml Lock	0001151-45	10

A series of horizontal dashed lines for writing notes.

RF Ablation


*Needles for RF and
Pulsed RF Ablation*



RFTL

Straight Thermolesion (RF) Needle with Active Tip

- 1 Straight Tip Geometry
- 2 Bevel Tip with Active Tip
- 3 Smooth Needle Coating



Product advantages at a glance:


- Reduced pain during puncture and needle placement
- Smooth gliding properties

Description	Item no.	PU	Description	Item no.	PU
22G x 52.6 mm (2") needle with 4mm Active Tip	001165-10B	25	20G x 52.6 mm (2") needle with 5mm Active Tip	001165-11D	25
22G x 98.6 mm (4") needle with 2.5mm Active Tip	001165-20A	25	20G x 98.6 mm (4") needle with 5mm Active Tip	001165-21B	25
22G x 98.6 mm (4") needle with 5mm Active Tip	001165-20B	25	20G x 98.6 mm (4") needle with 10mm Active Tip	001165-30D	25
22G x 98.6 mm (4") needle with 10mm Active Tip	001165-20C	25	20G x 142 mm (5 1/2") needle with 2.5mm Active Tip	001165-30A	10
			20G x 142 mm (5 1/2") needle with 5mm Active Tip	001165-30B	10
			20G x 142 mm (5 1/2") needle with 10mm Active Tip	001165-30C	10
			20G x 145 mm (5 3/4") needle with 5mm Active Tip	001165-30BM	10
			18G x 98.6 mm (4") needle with 10mm Active Tip	001165-22C	25

RFTL Curved

Curved Thermolesion (RF) Needle with Active Tip

- 1 Curved Tip Geometry
- 2 Bevel Tip with Active Tip
- 3 Smooth Needle Coating



Product advantages at a glance:


- Reduced pain during puncture and needle placement
- Smooth gliding properties

Description	Item no.	PU
22G x 98.6 mm (4") needle with 10mm Active Tip	011165-20C	25
20G x 98.6 mm (4") needle with 10mm Active Tip	011165-30D	25
18G x 98.6 mm (4") needle with 10mm Active Tip	011165-40D	25

RFTL Sono

Echogenic Thermolesion (RF) Needle with Active Tip

- 1 Cornerstone Reflectors
- 2 Bevel Tip with Active Tip
- 3 Smooth Needle Coating



Product advantages at a glance:


- Reliable and optimised needle echogenicity at higher angle
- Optimised ultrasound visibility of needle shaft
- Reduced pain during puncture and needle placement
- Smooth gliding properties

Description	Item no.	PU
22G x 98.6 mm (4") needle with 5mm Active Tip	001175-20B	25
20G x 52.6 mm (2") needle with 5mm Active Tip	001175-11D	25

PulsedRF

Thermolesion (RF) Needle with Active Tip

- 1 Bevel Tip with Active Tip
- 2 Hubless Design
- 3 Fixed Injection Tubing
- 4 Smooth Needle Coating



Product advantages at a glance:

- Reduced pain during puncture and needle placement
- Smooth gliding properties
- Reduced shadowing and improved visibility of target area
- Convenience due to closed system

Description	Item no.	PU
22G x 60 mm (2 1/4") needle with 5mm Active Tip, with Ø 1.6mm Plug	001255-12	10
22G x 100 mm (4") needle with 5mm Active Tip, with Ø 1.6mm Plug	001255-22	10
20G x 150 mm (6") needle with 5mm Active Tip, with Ø 1.6mm Plug	001255-32	10

AlliClip

Alligator Clip for PulsedRF Needles



Description	Item no.	PU
AlliClip Alligator Clip for PulsedRF Needles with \varnothing 1.6 mm Plug	001151-39D	10

A series of horizontal dashed lines for writing notes.

NOTES

A series of horizontal dashed lines for writing notes.

PAJUNK®

PAJUNK® GmbH Medizintechnologie
International Sales
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 9291-0
Fax: +49 7704 9291-600
info@pajunk.com
pajunk.com

PAJUNK® Medical Produkte GmbH
Central Europe
Pajunkstraße 2
78187 Geisingen, Germany
Tel.: +49 7704 8008-0
Fax: +49 7704 8008-150
info@pajunk-medipro.com
pajunk.de & pajunk.eu

PAJUNK® Swiss AG
Switzerland
Neuhofstrasse 5A
6340 Baar, Switzerland
Tel.: +41 41 760 96 27
info@pajunk.ch
pajunk.ch

PAJUNK® Medical System L.P.
United States of America
4575 Marconi Drive
Alpharetta, GA 30005, USA
Tel.: +1 770.493.6832
Fax: +1 678.514.3388
info@pajunk-usa.com
pajunkusa.com

PAJUNK® UK Medical Products Ltd
United Kingdom
Unit D1 The Waterfront
Goldcrest Way, Newburn Riverside
Newcastle upon Tyne, NE15 8NY, UK
Tel.: +44 191 264 7333
info@pajunk.co.uk
pajunk.co.uk